

M3 Vital Signs Monitor



Outpatient



Transport



Ward



Portable Design



Infrared Ear Thermometer



M3

Vital Signs Monitor

Proven to be a durable and efficient solution for outpatients and general wards, the M3 brings effective vital signs monitoring to most ambulatory and regular healthcares.

- 5.6" high resolution color TFT-LCD screen
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 8 hours of continuous work
- USB flash disk support

30,000

NIBP Measurement Review

72 h


Trend Review

800

Alarm Review

- Barcode scanner support
- Bi-directional communication with MFM-CMS central monitoring system
- Compatible with PatientCare Viewer PC data management software
- Wired LAN/Wi-Fi support, HL7 support via XML files or MFM-CMS
- Nurse call
- Suitable for adult, pediatric and neonatal patients

Seamless Connectivity



Mounting Solutions



Wall Mount



Rolling Stand



Bed Rail Clamp

Specification

Configuration	SpO ₂	NIBP	Quick Temp	Infrared Ear Thermometer
M3 (SpO ₂ only)	✓	—	—	—
M3 (NIBP only)	—	✓	—	—
M3 (NIBP+Quick Temp)	—	✓	✓	—
M3 (NIBP+Infrared Ear Thermometer)	—	✓	—	✓
M3 (SpO ₂ +NIBP+Quick Temp)	✓	✓	✓	—
M3 (SpO ₂ +NIBP+Infrared Ear Thermometer)	✓	✓	—	✓

4204 Jutland Drive,
Suite B, San Diego, CA 92117
Toll Free Phone Number:
888.850.4597