

# Scavenger System

## The Porter "Double Mask" Scavenger

### Unique, Comfortable, Quiet and Safe

The Porter scavenging system is the only double mask (mask-within-a-mask) system. It features an outer mask and a flexible, disposable (single-use), or autoclavable (multi-use), easily removable inner liner. It is the soft inner liner that creates the complete, yet comfortable, seal around the patient's nose, thus assuring the Porter mask's superior scavenging efficiency. The Porter masks are available in both adult and pedo sizes; the inner liners either unscented or citrus scented. The masks are autoclavable by either steam or chemical sterilizing and the Porter scavenging system is retrofittable to most flowmeters.



**Rated #1 in Scavenging Efficiency in Independent Clinical Trial.** In a major Dental School\* study of the effectiveness of different scavenging systems used in actual pediatric dental treatment, the Porter "Double Mask" Scavenging System was found to be the only one tested to reduce the occupational exposure of nitrous oxide to acceptable levels, meeting the 25 ppm requirements set by the National Institute of Occupational Safety and Health.

\* University of Medicine and Dentistry of New Jersey, New Jersey Dental School



Green Unscented Disposable



White Citrus Scent Disposable



Gray Autoclavable

Help promote office efficiency and infection control, these disposable or autoclavable latex-free liners can be on and gone, while complying fully with infection control protocols. Available in both adult and pedo sizes, unscented and citrus scented.

