

FEATURES

Part # GMF-LS01

Ports: ACMI, WOLF, STORZ, OLYMPUS

Dimensions: L: 12" W: 10" H: 6"

Weight: 16 lbs.

Power Input: 100/240 VAC 50/60 Hz

Circuit Protection: 8 AMP @ 250V

Safety Compliances: UL 2601-1

FDA 510K Listed

Lamp: 300 Watt Short ARC Xenon

Wattage: 300 W

Color Temperature: 6000° Kelvin

Light Source Warranty: 2 years

Lamp Warranty: 500 hours



Multi port turret, front panel includes light output LED and re-settable lamp timer.



OPTIONAL ACCESSORY

Part # GMF-LSST
Moveable stand



©Copyright 2005 Gulf Medical Fiberoptics, Inc.