


PEDIGO
fine healthcare equipment

[Home](#)
[Products](#)
[Stretchers](#)
[About Pedigo](#)
[Find a Dealer](#)
[Find a Rep](#)
[Show Schedule](#)
[Contact Us](#)

Footstools and Step Stools

- ▶ Strong, durable, all welded construction
- ▶ Legs and handrails constructed of heavy 16-gauge 1" tube
- ▶ Slip resistant corrugated rubber top surface
- ▶ Rubber floor tips
- ▶ Tip resistant design

[P-10](#)
[P-1011-SS](#)
[P-1010-A-SS](#)
[P-11](#)
[P-1012-SS](#)
[P-18](#)
[P-1010-SS](#)
[P-10-A](#)
[P-18-A](#)


Footstools

Chrome Models

P-10

- ▶ Dimensions: 12"D x 16"W x 8-1/4"H

P-11

- ▶ Dimensions: 12"D x 30"W x 8"H

Stainless Steel Models

P-1010-SS

- ▶ Dimensions: 12"D x 16"W x 8-1/4"H

P-1011-SS

- ▶ Dimensions: 12"D x 30"W x 5"H

P-1012-SS

- ▶ Dimensions: 12"D x 30"W x 8"H



Footstools with Handrail

Chrome Models

P-10-A

- ▶ Same as P-10 except with handrail
- ▶ Curved handle extends 38-1/2" above the floor

Stainless Steel Models

P-1010-A-SS

- ▶ Same as P-1010-SS except with handrail
- ▶ Curved handle extends 38-1/2" above the floor



Double Step Stools

Chrome Models

P-18

- ▶ Steps are 12"D x 16"W and 8" and 16" high

P-18-A

- ▶ Same as P-18 except with handrail
- ▶ Handrail is 47" above the floor

[Home](#) - [Top](#) - [Products](#) - [About Pedigo](#) - [Contact Us](#)

Pedigo Products, Inc.
4000 SE Columbia Way
Vancouver WA, 98661 USA

Phone: (800) 822-3501
Fax: (360) 696-1700

Copyright © 2000 Pedigo Products, Inc. All Rights Reserved.