


**PEDIGO**
*fine healthcare equipment*

[Home](#)
[Products](#)
[Stretchers](#)
[About Pedigo](#)
[Find a Dealer](#)
[Find a Rep](#)
[Show Schedule](#)
[Contact Us](#)

## Mayo Tables

- ▶ All welded construction for strength and durability
- ▶ Choice of stainless steel, chrome plated or our new, durable *Near-Chrome* powder-coated finish models
- ▶ 2" double ball bearing swivel casters

[P-1065-SS](#)
[P-1068-SS](#)
[P-1069-A-SS](#)
[P-1066-SS](#)
[P-1068-A-SS](#)
[P-66](#)
[P-1066-A-SS](#)
[P-1069-SS](#)
[P-65](#)


### **P-1065-SS**

- ▶ All welded stainless steel construction
- ▶ Removable, 20" x 25" stainless steel tray for easy cleaning
- ▶ Foot-pedal operated friction lock adjusts height from 36" to 62"
- ▶ Base designed to move under low clearance equipment

### **P-1066-SS**

- ▶ Same as P-1065-SS except for tray dimensions
- ▶ 16-1/4" x 21-1/4" tray

### **P-1066-A-SS**

- ▶ Same as P-1066-SS except with hand-operated height adjustment

### **P-1068-SS**

- ▶ Same as P-1065-SS except for tray dimensions
- ▶ 12-5/8" x 19-1/8" tray

### **P-1068-A-SS**

- ▶ Same as P-1068-SS except with hand-operated height adjustment



### **P-1069-SS**

- ▶ All welded stainless steel construction
- ▶ Removable, 16-1/4" x 21-1/4" stainless steel tray
- ▶ One forward leg with caster for increased mobility
- ▶ Foot-pedal operated friction lock adjusts height from 36" to 62"
- ▶ Base designed to move under low clearance equipment

### **P-1069-A-SS**

- ▶ Same as P-1069-SS except with hand-operated height adjustment

### **P-66**

- ▶ Removable, 19-1/8" x 12-5/8" stainless steel tray on *Near-Chrome* powder-coated frame
- ▶ Friction lock adjusts height from 32" to 53"

**P-65**

- ▶ Removable, 19-1/8" x 12-5/8" stainless steel tray on chrome plated frame
- ▶ Friction lock adjusts height from 30" to 40"

[Home](#) - [Top](#) - [Products](#) - [About Pedigo](#) - [Contact Us](#)

---

Pedigo Products, Inc.  
4000 SE Columbia Way  
Vancouver WA, 98661 USA

**Phone:** (800) 822-3501  
**Fax:** (360) 696-1700

Copyright © 2000 Pedigo Products, Inc. All Rights Reserved.