



PEDIGO

fine healthcare equipment



Home

Products

Stretchers

About Pedigo

Find a Dealer

Find a Rep

Show Schedule

Contact Us

Small Utility Tables

- ▶ Easy to Clean
- ▶ Stainless steel, Steril-Gard construction eliminates full sub-tops
- ▶ Improves the control of undesirable substances
- ▶ Sound-proofed reinforcing channels for quiet use

▶ Available Models:

[SG-80-A-SS](#)

[SG-80-SS](#)

[SG-81-SS](#)



SG-80-A-SS

- ▶ Utility table with shelf, rail and two drawers
- ▶ 1-1/2" H, 3-sided top rail
- ▶ 3" double ball bearing, rubber swivel casters for silent mobility
- ▶ Dimensions: 16"W x 20"L x 34"H; Drawers 12-1/2"L x 15-1/2"W x 3-1/2"D

SG-80-SS

- ▶ Utility table with shelf, rail and one drawer
- ▶ 1-1/2" H, 3-sided top rail
3" double ball bearing, rubber swivel casters for silent mobility
- ▶ Dimensions: 16"W x 20"L x 34"H; Drawer 12-1/2"L x 15-1/2"W x 3-1/2"D

**SG-81-SS**

- ▶ Utility table with shelf and rail, no drawer
- ▶ 1-1/2" H, 3-sided top rail
- ▶ 3" double ball bearing, rubber swivel casters for silent mobility
- ▶ Dimensions: 16"W x 20"L x 34"H

[Home](#) - [Top](#) - [Products](#) - [About Pedigo](#) - [Contact Us](#)

Pedigo Products, Inc.
4000 SE Columbia Way
Vancouver WA, 98661 USA

Phone: (800) 822-3501
Fax: (360) 696-1700

Copyright © 2000 Pedigo Products, Inc. All Rights Reserved.