

# tuttnauer

Sterilization & Infection Control Products

## *Closed Door System Fully Automatic Autoclaves*



New **EX9™** and **EX10™**  
autoclaves will meet all your sterilization needs

- Chamber door remains closed during entire cycle
- Hepa filtered air is circulated in chamber during drying
- All program parameters can be customized and stored
- New brighter and easier to read panel
- Two year parts and labor warranty

\* Above photo shown with optional printer

**W**ith the simplicity of one touch design all your sterilization and drying needs are fulfilled. The Tuttnauer fully automatic autoclaves fill, sterilize, exhaust and dry at the touch of a button.

## Closed Door System Fully Automatic Autoclaves

### EZ9™ and EZ10™ Series

Item No.	Chamber Dimensions Ø x L In(mm)	Volume	Overall Dimensions DxWxH In(mm)	Standard Cassettes	Tray Dimensions DxWxH In(mm)	No. of Trays	Standard Unwrapped Cycle time	Voltage(V) Freq.(Hz)	Power(W) Current(A)	Shipping Weight lbs(kg)
EZ9	9" x 18" (230 x 470)	5gal/19L	21.5"x20"x14.4" (545x510x365)	2 full/2 half	6.7"x16.3"x0.8" (170x415x20)	3	Cold: 30 mins Hot: 20 mins	120V 50/60	1400W 12A	95 (43)
EZ10	10" x 19" (254 x 475)	6gal/23L	21.5"x20"x14.4" (545x510x365)	3 full/3 half	6.7"x16.3"x0.8" (170x415x20)	4	Cold: 30 mins Hot: 20 mins	120V 50/60	1400W 12A	95 (43)
EZ10k Kwiklave	10" x 19" (254 x 475)	6gal/23L	21.5"x20"x14.4" (545x510x365)	3 full/3 half	6.7"x16.3"x0.8" (170x415x20)	4	Cold: 12 mins Hot: 9 mins	230V 50/60	2200W 10A	95 (43)

Model EZ10k requires a voltage between 220V and 235V for proper operation. When applicable a buck/boost transformer (0.5KVA) may be required.

#### OPTIONAL INTERNAL PRINTER



documents:  
date, time,  
temperature,  
pressure.

#### CONTROL PANEL



1. Program selection
2. Detailed digital display
3. Programmable keys
4. System indicators

New brighter and easier to read panel


#### FEATURES:

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches.
- High efficiency air pump circulates air in chamber during drying cycle.
- 0.2µm Hepa air filter (replaceable)
- New brighter and easier to read panel.
- Program parameters can be customized and stored.

#### Clean & Simple Ultrasonic Cleaners and Tablets



#### Instrument Pouch Rack



#### Steam Distiller



#### Chamber Brite Autoclave Cleaner



[www.tuttnauer.com](http://www.tuttnauer.com)

**tuttnauer**  
Sterilization & Infection Control Products

#### Tuttnauer USA Co.

25 Power Drive Hauppauge NY, 11788, USA  
Tel: (800) 624 5836, (631) 737 4850 Fax: (631) 737 0720  
Email: [info@tuttnauer.com](mailto:info@tuttnauer.com)