

**NEW!**

**smiths medical**  
bringing technology to life

## Medfusion® 3500 Syringe Pump with PharmGuard® v5 Software

Medfusion® 3500  
syringe pump  
now recognizes

**Enteral  
Syringes**



**For safe, accurate, enteral  
nutrition delivery**

**medfusion**

™

# Reduce the risk of IV / enteral misconnections

## Medfusion® 3500 syringe pumps with PharmGuard® v5 software now recognize NeoMed®, Baxa® NeoThrive™, and Vygon® enteral syringes

Adding enteral syringe data to PharmGuard® v5 software allows you to:

- Meet Joint Commission requirements for enteral feedings<sup>1</sup>
- Differentiate your patients' IV and enteral lines<sup>2</sup>
- Meet your facility's needs and provide smarter, safer enteral delivery

### Supported Enteral Syringes

NeoMed® - 2 sizes - 20 mL, 60 mL - orange  
Baxa® NeoThrive™ - 3 sizes - 20 mL, 35 mL, 60 mL - orange  
Vygon® C-Gon - 2 sizes - 20 mL, 60 mL - purple

## Smart and safe just got smarter and safer

By partnering with the leading enteral syringe companies in the market, Smiths Medical ensures you'll have the best of both worlds - a smarter infusion system and safer enteral delivery

- Orange or purple coloring for easy enteral identification
- Unique tips will not connect to standard I.V. Luers



**neothrive**  
Brought to you by Baxa Corporation.

**VYCON**

**NEOMED**  
incorporated

Call Smiths Medical today to find out how to incorporate enteral syringe data into your Medfusion® 3500 pumps with PharmGuard® v5 software

**1-800-258-5361**

## Order Information

PharmGuard® Supported Syringes Series 1 67-2449-51

1. Guenter, Peggi, PhD., R.N., C.N.S.N., et al. *Enteral Feeding Misconnections: A Consortium Position Statement*. The Joint Commission Journal on Quality and Patient Safety. May 2008; Volume 34 Number 5
2. Abstract: *Enteral feeding set connectors and adapters*. Association for the Advancement of Medical Instrumentation. ANSI/AAMI ID54:1996/©2005.

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

### Smiths Medical ASD, Inc.

St. Paul, MN 55112, USA  
Phone: 1-214-618-0218  
Toll-Free USA: 1-800-258-5361

www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc.

Product(s) described may not be licensed or available for sale in Canada or other countries outside of the United States.

Medfusion, PharmGuard and the Smiths Medical and Medfusion design marks are trademarks of the Smiths Medical family of companies. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners.  
©2009 Smiths Medical family of companies. All rights reserved. 10/09 193816

### Smiths Medical Canada Ltd.

Markham, Ontario, Canada, L3R 4Y8  
Phone: 905-477-2000  
Toll-Free: 1-800-387-4346

### Customer and Clinical Services

U.S. 1-800-258-5361  
Canada 1-800-387-4346

www.smiths-medical.com

**smiths medical**